


Anti-Src Rabbit Polyclonal Antibody For Wb And Elisa - P12931

Item Number: CM0015708



Introduction

This rabbit polyclonal antibody targets human Src (c-Src, p60-Src) with cross-reactivity to mouse and rat. Validated for WB and ELISA, it supports research into tyrosine kinase signaling, cell adhesion, and oncogenesis.

[Learn More](#)

Parameter	Value
Product Name	Src Rabbit pAb
Remarks/Alias	ASV; SRC1; THC6; c-SRC; p60-Src; Src
Species	Human
Gene ID (Human)	6714
Gene ID	6714
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 410-536 of human Src (NP_005408.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P12931
Protein Weight	60kDa
Shipping	Ice bag