

Anti-Srd5A1 Polyclonal Antibody For Wb, Elisa - P18405

Item Number: CM0015001

Introduction

Rabbit polyclonal antibody against human SRD5A1, validated for WB and ELISA. Cross-reacts with mouse and rat. Ideal for androgen signaling, steroid metabolism, and sexual differentiation research.

[Learn More](#)

Parameter	Value
Product Name	SRD5A1 Rabbit pAb
Remarks/Alias	S5AR 1; SRD5A1
Species	Human
Gene ID	6715
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 100-200 of human SRD5A1 (NP_001038.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P18405
Protein Weight	29kDa
Shipping	Ice bag