


Anti-Srd5A2 Monoclonal Antibody For Wb, Elisa - P31213

Item Number: CM0003682



Introduction

High-quality rabbit monoclonal antibody targeting SRD5A2, validated for WB and ELISA. Detects human, mouse, and rat SRD5A2, a key enzyme in androgen metabolism. Suitable for prostate biology research. P31213.

[Learn More](#)

Parameter	Value
Product Name	SRD5A2 Rabbit mAb
Species	Human
GeneID	6716 (Human)
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human SRD5A2 (P31213).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P31213
Protein Weight	28 kDa
Shipping	Ice bag