

Anti-Srebf1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P36956

Item Number: CM0016375



Introduction

High-quality rabbit polyclonal antibody targeting human SREBF1 (SREBP1), validated for WB, IF/ICC, and ELISA. Ideal for studying cholesterol biosynthesis, lipid homeostasis, and metabolic regulation. Cross-reacts with mouse and rat.

[Learn More](#)

Parameter	Value
Product Name	SREBF1 Rabbit pAb
Remarks/Alias	HMD, IFAP2, SREBP1, bHLHd1, SREBF1
Species	Human
GeneID (Human)	6720
GeneID	6720
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 200-310 of human SREBF1 (NP_004167.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P36956
Protein Weight	122kDa
Shipping	Ice bag