


Anti-Srebf2 Rabbit Polyclonal Antibody For Wb, If/Icc, If-P, Elisa - Q12772

Item Number: CM0020421



Introduction

Rabbit polyclonal antibody targeting human SREBF2 (SREBP-2), a master cholesterol biosynthesis transcription factor. Suitable for Western blot, immunofluorescence, and ELISA; cross-reacts with mouse. For lipid metabolism studies.

[Learn More](#)

Parameter	Value
Product Name	SREBF2 Rabbit pAb
Remarks/Alias	SREBP2; bHLHd2; SREBP-2
Species	Human
Gene ID (Human)	6721
Gene ID	6721
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 401-500 of human SREBF2 (NP_004590.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IF-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q12772
Protein Weight	124kDa
Shipping	Ice bag