


Anti-Srebp2 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q12772

Item Number: CM0015831



Introduction

SREBP2 Rabbit Polyclonal Antibody targeting human sterol regulatory element-binding protein 2 (Q12772). Suitable for WB, IHC-P, ELISA; cross-reacts with mouse and rat. Essential tool for studying cholesterol biosynthesis and lipid metabolism pathways.

[Learn More](#)

Parameter	Value
Product Name	SREBP2 Rabbit pAb
Remarks/Alias	SREBP2; bHLHd2; SREBP-2
Species	Human
GeneID (Human)	6721
GeneID	6721
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 155-310 of human SREBF2 (NP_004590.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q12772
Protein Weight	124kDa
Shipping	Ice bag