

Anti-Srf Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P11831

Item Number: CM0026216



Introduction

High-quality rabbit polyclonal antibody against human SRF, validated for WB, IHC-P, ELISA with cross-reactivity to mouse and rat. Suitable for transcription regulation and cardiac development research.

[Learn More](#)

Parameter	Value
Product Name	SRF Rabbit pAb
Remarks/Alias	MCM1; SRF
Species	Human
Gene ID (Human)	6722
Gene ID	6722
Immunogen	Synthetic peptide; A synthetic peptide corresponding to a sequence within amino acids 321-420 of human SRF (NP_003122.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P11831
Protein Weight	52kDa
Shipping	Ice bag