

Anti-Srfbp1 Rabbit Polyclonal Antibody For Wb, Elisa - Q8Nef9

Item Number: CM0027341

Introduction

Rabbit polyclonal antibody against human SRFBP1, validated for Western blot (WB) and ELISA, cross-reactive with mouse. Ideal for studying cardiac gene transcription and SLC2A4 processing.

[Learn More](#)

Parameter	Value
Product Name	SRFBP1 Rabbit pAb
Remarks/Alias	P49; Rib1; BUD22; STRAP; p49/STRAP; SRFBP1
Species	Human
Gene ID	153443
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human SRFBP1 (NP_689759.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q8NEF9
Protein Weight	49kDa
Shipping	Ice bag