

Anti-Srgn Polyclonal Antibody For Wb, Elisa - P10124

Item Number: CM0013596

Introduction

High-quality rabbit polyclonal antibody against human SRGN, validated for WB and ELISA. Detects both human and mouse serglycin. Suitable for research on secretory granules, immune cell function, and bone mineralization inhibition.

[Learn More](#)

Parameter	Value
Product Name	SRGN Rabbit pAb
Remarks/Alias	PPG; PRG; PRG1; SRGN
Species	Human
GeneID	5552
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 28-158 of human SRGN (NP_002718.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P10124
Protein Weight	18kDa
Shipping	Ice bag