

## Anti-Srgn Polyclonal Antibody For Wb, If-P, Elisa - P10124

Item Number: CM0028082

### Introduction

High-quality rabbit polyclonal antibody targeting human SRGN (Serglycin). Suitable for WB, IF-P, ELISA; cross-reacts with human and mouse. Ideal for apoptosis and secretory granule research.

[Learn More](#)

| Parameter        | Value                                                                             |
|------------------|-----------------------------------------------------------------------------------|
| Product Name     | SRGN Rabbit pAb                                                                   |
| Remarks/Alias    | PPG; PRG; PRG1                                                                    |
| Species          | Human                                                                             |
| GeneID (Human)   | 5552                                                                              |
| GeneID           | 5552                                                                              |
| Immunogen        | Synthetic peptide corresponding to amino acids 59-158 of human SRGN (NP_002718.2) |
| Source           | Rabbit                                                                            |
| Category         | Polyclonal Antibodies                                                             |
| Application      | WB, IF-P, ELISA                                                                   |
| Cross Reactivity | Human, Mouse                                                                      |
| SWISS            | P10124                                                                            |
| Protein Weight   | 18kDa                                                                             |
| Shipping         | Ice bag                                                                           |