

Anti-Sri Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P30626

Item Number: CM0030718



Introduction

KO-validated rabbit polyclonal antibody against SRI (sorcin), a cardiac calcium-binding protein. Suitable for WB, IHC-P, and ELISA with cross-reactivity to human, mouse, and rat. Ideal for excitation-contraction coupling research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] SRI Rabbit pAb
Remarks/Alias	SCN; V19; CP22; CP-22; RI
Species	Human
GeneID (human)	6717
GeneID	6717
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-198 of human SRI (NP_003121.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P30626
Protein Weight	22kDa
Shipping	Ice bag