


Anti-Srm Polyclonal Antibody For Wb, Elisa - P19623

Item Number: CM0029797



Introduction

SRM Rabbit pAb is a polyclonal antibody targeting spermidine synthase, validated for WB and ELISA with cross-reactivity to human, mouse, and rat. Suitable for polyamine biosynthesis studies.

[Learn More](#)

Parameter	Value
Product Name	SRM Rabbit pAb
Remarks/Alias	PAPT; SPS1; SPDSY; SRML1; SRM
Species	Human
GeneID (Human)	6723
GeneID	6723
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-302 of human SRM (NP_003123.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P19623
Protein Weight	34kDa
Shipping	Ice bag