


Anti-Srp14 Polyclonal Antibody For Wb, Elisa - P37108

Item Number: CM0015595



Introduction

Rabbit polyclonal antibody targeting human SRP14, validated for WB and ELISA. Recognizes human, mouse, and rat SRP14, a 15 kDa signal recognition particle component. Ideal for studying protein translocation to the ER.

[Learn More](#)

Parameter	Value
Product Name	SRP14 Rabbit pAb
Remarks/Alias	ALURBP; SRP14
Species	Human
GeneID	6727
Immunogen	Recombinant fusion protein corresponding to amino acids 1-136 of human SRP14 (NP_003125.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P37108
Protein Weight	15kDa
Shipping	Ice bag