


Anti-Srp19 Polyclonal Antibody For Wb, Ihc-P, Elisa - P09132

Item Number: CM0030774



Introduction

Rabbit polyclonal antibody targeting SRP19, a component of the signal recognition particle. Validated for WB, IHC-P, and ELISA with reactivity across human, mouse, and rat samples. Ideal for protein trafficking research.

[Learn More](#)

Parameter	Value
Product Name	SRP19 Rabbit pAb
Remarks/Alias	SRP19
Species	Human
Gene ID	6728
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-144 of human SRP19 (NP_003126.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P09132
Protein Weight	16kDa
Shipping	Ice bag