


Anti-Srp54 Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - P61011

Item Number: CM0006951



Introduction

Rabbit monoclonal antibody against human SRP54 (P61011), validated for WB, IF/ICC, ELISA. Reacts with human, mouse, rat. SRP54 is a signal recognition particle subunit involved in cotranslational protein targeting to ER.

[Learn More](#)

Parameter	Value
Product Name	SRP54 Rabbit mAb
Remarks/Alias	SCN8; SRP54
Species	Human
Gene ID (Human)	6729
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 400-504 of human SRP54 (P61011).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P61011
Protein Weight	56kDa
Shipping	Ice bag