


Anti-Srp72 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - O76094

Item Number: CM0017139



Introduction

High-quality rabbit polyclonal antibody against human SRP72, validated for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with mouse and rat. Ideal for studying SRP-mediated protein targeting to the ER.

[Learn More](#)

Parameter	Value
Product Name	SRP72 Rabbit pAb
Remarks/Alias	BMFF; BMFS1; HEL103; SRP72
Species	Human
Gene ID (Human)	6731
Gene ID	6731
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 440-544 of human SRP72 (NP_008878.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O76094
Protein Weight	75kDa
Shipping	Ice bag