

Anti-Srp72 Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - O76094

Item Number: CM0012720



Introduction

High-quality rabbit monoclonal antibody against human SRP72 (O76094). Validated for WB, IHC-P, IF/ICC, ELISA. Reacts with human, mouse, rat. Ideal for studying SRP complex and ER targeting.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] SRP72 Rabbit mAb
Remarks/Alias	BMFF; BMFS1; HEL103; SRP72
Species	Human
GeneID (Human)	6731
GeneID	6731
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 440-544 of human SRP72 (NP_008878.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O76094
Protein Weight	75kDa
Shipping	Ice bag