


Anti-Srp9 Polyclonal Antibody For Wb, If/Icc, Elisa - P49458

Item Number: CM0030997



Introduction

Anti-SRP9 Rabbit Polyclonal Antibody validated for WB, IF/ICC, ELISA. Recognizes human and mouse SRP9, a signal recognition particle component mediating ER targeting. Recombinant immunogen (aa 1-86).

[Learn More](#)

Parameter	Value
Product Name	SRP9 Rabbit pAb
Remarks/Alias	ALURBP; SRP9
Species	Human
GeneID (Human)	6726
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-86 of human SRP9 (NP_003124.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P49458
Protein Weight	10kDa
Shipping	Ice bag