


Anti-Srpk1 Monoclonal Antibody For Wb, Ihc-P, Elisa - Q96Sb4

Item Number: CM0010514



Introduction

High-quality rabbit monoclonal anti-SRPK1 antibody, validated for Western blot, IHC-P, and ELISA. Reacts with human, mouse, rat SRPK1 (UniProt Q96SB4). Ideal for splicing and cell cycle research.

[Learn More](#)

Parameter	Value
Product Name	SRPK1 Rabbit mAb
Remarks/Alias	SFRSK1
Species	Human
Gene ID	6732
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 291-390 of human SRPK1 (NP_003128.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96SB4
Protein Weight	74kDa/72kDa/92kDa
Shipping	Ice bag