

Anti-SrpK1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q96Sb4

Item Number: CM0025583



Introduction

SRPK1 rabbit polyclonal antibody (UniProt Q96SB4) validated for Western blot, immunofluorescence, and ELISA. Cross-reacts with human, mouse, rat. Suitable for studying SR protein phosphorylation, splicing regulation, and cell cycle progression.

[Learn More](#)

Parameter	Value
Product Name	SRPK1 Rabbit pAb
Remarks/Alias	SFRSK1; SRPK1
Species	Human
GeneID (Human)	6732
GeneID	6732
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 290-480 of human SRPK1 (NP_003128.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96SB4
Protein Weight	74kDa
Shipping	Ice bag