

Anti-Srpr Rabbit Polyclonal Antibody For Wb, Elisa - P08240

Item Number: CM0027180



Introduction

Rabbit polyclonal antibody targeting human SRPR (SRP receptor alpha subunit) for Western blot and ELISA applications. Cross-reacts with human; immunogen covers N-terminal region. For research use.

[Learn More](#)

Parameter	Value
Product Name	SRPR Rabbit pAb
Remarks/Alias	DP; SRPR; Sralpha
Species	Human
GeneID (Human)	6734
GeneID	6734
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P08240
Protein Weight	70kDa
Shipping	Ice bag