


Anti-Srprb Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Y5M8

Item Number: CM0015875



Introduction

SRPRB Rabbit Polyclonal Antibody validated for WB, IF/ICC, ELISA in human, mouse, rat. Targets signal recognition particle receptor subunit beta for protein targeting studies.

[Learn More](#)

Parameter	Value
Product Name	SRPRB Rabbit pAb
Remarks/Alias	APMCF1; SR-beta; SRPRB
Species	Human
Gene ID (Human)	58477
Gene ID	58477
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 57-271 of human SRPRB (NP_067026.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y5M8
Protein Weight	30kDa
Shipping	Ice bag