

Anti-Srpx Polyclonal Antibody For Wb, Ihc-P, Elisa - P78539

Item Number: CM0014832

Introduction

Human SRPX rabbit polyclonal antibody validated for Western blot, IHC-P, and ELISA. Cross-reacts with mouse. Ideal for studying phagocytosis, cell adhesion, and retinal disk shedding.

[Learn More](#)

Parameter	Value
Product Name	SRPX Rabbit pAb
Remarks/Alias	DRS; ETX1; SRPX1; HEL-S-83p; SRPX
Species	Human
GeneID	8406
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 40-320 of human SRPX (NP_006298.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P78539
Protein Weight	52kDa
Shipping	Ice bag