


Anti-Srrm1 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q8Iyb3

Item Number: CM0029939



Introduction

Rabbit polyclonal antibody against human SRRM1 (SRm160), suitable for WB, IHC-P, ELISA. Cross-reacts with mouse. Ideal for RNA splicing and spliceosome research.

[Learn More](#)

Parameter	Value
Product Name	SRRM1 Rabbit pAb
Remarks/Alias	160-KD; POP101; SRM160; SRRM1
Species	Human
GeneID (human)	10250
GeneID	10250
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-160 of human SRRM1 (NP_005830.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q8IYB3
Protein Weight	102kDa
Shipping	Ice bag