

Anti-Srrt Rabbit Polyclonal Antibody For Ihc-P, Elisa - Q9Bxp5

Item Number: CM0028749



Introduction

High-quality SRRT/ARS2 rabbit polyclonal antibody for IHC-P and ELISA. Detects human, mouse, rat SRRT protein (101 kDa, Q9BXP5). Essential tool for investigating miRNA processing, cell cycle progression, and transcription regulation.

[Learn More](#)

Parameter	Value
Product Name	SRRT Rabbit pAb
Remarks/Alias	ARS2, ASR2, serrate, SRRT
Species	Human
GeneID (Human)	51593
GeneID	51593
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 623-872 of human SRRT (NP_001122325.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BXP5
Protein Weight	101kDa
Shipping	Ice bag