


Anti-Srsf1/Sf2/Asf Monoclonal Antibody For Wb, Ihc-P, Elisa - Q07955

Item Number: CM0009576



Introduction

High-quality rabbit monoclonal antibody targeting SRSF1 (SF2/ASF), validated for WB, IHC-P, and ELISA. Recognizes human, mouse, and rat SRSF1. Ideal for splicing research, diagnostic development, and biomarker studies.

[Learn More](#)

Parameter	Value
Product Name	SRSF1/SF2/ASF Rabbit mAb
Remarks/Alias	ASF; SF2; SFRS1; SF2p33; SRp30a; SRSF1/SF2/ASF
Species	Human
Gene ID (Human)	6426
Gene ID	6426
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 100-200 of human SRSF1/SF2/ASF (NP_008855.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q07955
Protein Weight	28kDa
Shipping	Ice bag