


Anti-Srsf2 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q01130

Item Number: CM0029531



Introduction

Rabbit polyclonal antibody targeting human SRSF2 splicing factor for Western blot (WB), IHC-P, and ELISA. Cross-reacts with human, mouse, and rat. Essential for pre-mRNA splicing research and diagnostic applications.

[Learn More](#)

Parameter	Value
Product Name	SRSF2 Rabbit pAb
Remarks/Alias	SC35; PR264; SC-35; SFRS2; SFRS2A; SRp30b; SRSF2
Species	Human
Gene ID (Human)	6427
Gene ID	6427
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human SRSF2 (NP_001182356.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q01130
Protein Weight	25kDa
Shipping	Ice bag