


Anti-Srsf3 Monoclonal Antibody For Wb, Ihc-P, Elisa - P84103

Item Number: CM0006265



Introduction

Rabbit recombinant monoclonal antibody targeting SRSF3, a serine/arginine-rich splicing factor involved in mRNA alternative splicing and nuclear export. Validated for WB, IHC-P, ELISA in human, mouse, and rat samples.

[Learn More](#)

Parameter	Value
Product Name	SRSF3 Rabbit mAb
Remarks/Alias	SFRS3; SRp20; SRSF3
Species	Human
Gene ID (Human)	6428
Gene ID	6428
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human SRSF3 (P84103)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P84103
Protein Weight	19kDa
Shipping	Ice bag