


Anti-Srsf3 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P84103

Item Number: CM0030262



Introduction

Rabbit polyclonal antibody targeting SRSF3 (SRp20), validated for WB, IHC-P, IF/ICC, and ELISA. Recognizes human, mouse, and rat. Suitable for splicing and mRNA export research.

[Learn More](#)

Parameter	Value
Product Name	SRSF3 Rabbit pAb
Remarks/Alias	SFRS3; SRp20; SRSF3
Species	Human
GeneID (Human)	6428
GeneID	6428
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-90 of human SRSF3 (NP_003008.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P84103
Protein Weight	19kDa
Shipping	Ice bag