

Anti-Srsf4 Polyclonal Antibody For Wb, Elisa - Q08170

Item Number: CM0029376

Introduction

Rabbit polyclonal antibody against human SRSF4 (SFRS4/SRP75, 57 kDa), validated for WB and ELISA with cross-reactivity to mouse and rat. Ideal for splicing research.

[Learn More](#)

Parameter	Value
Product Name	SRSF4 Rabbit pAb
Remarks/Alias	SFRS4; SRP75; SRSF4
Species	Human
GeneID (Human)	6429
GeneID	6429
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-180 of human SRSF4 (NP_005617.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q08170
Protein Weight	57kDa
Shipping	Ice bag