

Anti-Srsf7 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q16629

Item Number: CM0018804



Introduction

High-purity rabbit polyclonal antibody targeting human SRSF7/9G8, a splicing factor involved in mRNA processing and export. Validated for WB, IHC-P, IF/ICC, ELISA, with cross-reactivity to mouse and rat. Ideal for splicing regulation research.

[Learn More](#)

Parameter	Value
Product Name	SRSF7 Rabbit pAb
Remarks/Alias	9G8; AAG3; SFRS7; SRSF7
Species	Human
GeneID (human)	6432
GeneID	6432
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 60-120 of human SRSF7 (NP_001026854).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16629
Protein Weight	27kDa
Shipping	Ice bag