

Anti-Srsf8 Rabbit Polyclonal Antibody For Wb - Q9BrI6

Item Number: CM0021747



Introduction

Rabbit polyclonal anti-SRSF8 antibody validated for Western blot and ELISA, with mouse cross-reactivity. Targets human serine/arginine-rich splicing factor 8 involved in pre-mRNA alternative splicing.

[Learn More](#)

Parameter	Value
Product Name	SRSF8 Rabbit pAb
Remarks/Alias	DSM-1; SRP46; SFRS2B; SRSF8
Species	Human
Gene ID (Human)	10929
Gene ID	10929
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q9BRL6
Protein Weight	32kDa
Shipping	Ice bag