

Anti-Ssb Rabbit Polyclonal Antibody (Pab) For Wb, Ihc-P, If/Icc, Ip, Elisa - P05455

Item Number: CM0013018



Introduction

Validate SSB/La protein expression in human samples with this rabbit polyclonal antibody. Suitable for WB, IHC-P, IF/ICC, IP, and ELISA. Ideal for autoimmunity and RNA biology research.

[Learn More](#)

Parameter	Value
Product Name	SSB Rabbit pAb
Remarks/Alias	La; LARP3; La/SSB; SSB
Species	Human
GeneID (Human)	6741
GeneID	6741
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 166-408 of human SSB (NP_003133.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human
SWISS	P05455
Protein Weight	47kDa
Shipping	Ice bag