

Anti-Ssca1/Serpinb3 Monoclonal Antibody For Wb, If/Icc, Elisa - P29508

Item Number: CM0007760



Introduction

Monoclonal antibody targeting human Serpin B3/SCCA-1, validated for WB, IF/ICC, and ELISA. Cross-reacts with human, mouse, and rat. Ideal for cancer and apoptosis research.

[Learn More](#)

Parameter	Value
Product Name	SSCA1/SerpinB3 Rabbit mAb
Remarks/Alias	SCC; T4-A; SCCA1; SSCA1; SCCA-1; HsT1196; SCCA-PD; SSCA1/SerpinB3
Species	Human
Gene ID	6317 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 2-390 of human SSCA1/SerpinB3 (NP_008850.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P29508
Protein Weight	45kDa
Shipping	Ice bag