

Anti-Ssna1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - O43805

Item Number: CM0013671

Introduction

Rabbit polyclonal anti-SSNA1 antibody validated for WB, IHC-P, ELISA; reacts with human, mouse, rat. Targets microtubule nucleation factor involved in axon branching, cell division, and microtubule stability.

[Learn More](#)

Parameter	Value
Product Name	SSNA1 Rabbit pAb
Remarks/Alias	N14; NA14; NA-14; SSNA1
Species	Human
GeneID (human)	8636
GeneID	8636
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-119 of human SSNA1 (NP_003722.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O43805
Protein Weight	14kDa
Shipping	Ice bag