


Anti-Ssr3 Rabbit Polyclonal Antibody For Wb And Elisa - Q9Unl2

Item Number: CM0020618



Introduction

High-quality anti-SSR3 rabbit polyclonal antibody for Western blotting and ELISA. Cross-reacts with human, mouse, and rat SSR3. Reliable tool for studying ER protein retention and calcium-binding mechanisms.

[Learn More](#)

Parameter	Value
Product Name	SSR3 Rabbit pAb
Remarks/Alias	TRAPG; SSR3
Species	Human
Gene ID	6747
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 70-145 of human SSR3 (NP_009038.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UNL2
Protein Weight	21kDa
Shipping	Ice bag