


Anti-Ssrp1 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q08945

Item Number: CM0025291



Introduction

Rabbit anti-SSRP1 polyclonal antibody (pAb) validated for WB, IHC-P, IF/ICC, ELISA. Reacts with human and mouse SSRP1 (FACT80), a FACT complex subunit involved in chromatin remodeling, transcription elongation, and DNA repair.

[Learn More](#)

Parameter	Specification
Product Name	SSRP1 Rabbit pAb
Remarks/Alias	FACT; T160; FACT80; SSRP1
Species	Human
Gene ID (Human)	6749
Gene ID	6749
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human SSRP1 (NP_003137.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q08945
Protein Weight	81kDa
Shipping	Ice bag