


Anti-Ssrp1 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q08945

Item Number: CM0015178



Introduction

Rabbit polyclonal antibody against human SSRP1 (FACT80) validated in WB, IHC-P, IF/ICC, ELISA. Targets chromatin transcription elongation factor subunit; ideal for nucleosome remodeling, DNA repair, and transcriptional regulation research.

[Learn More](#)

Parameter	Value
Product Name	SSRP1 Rabbit pAb
Remarks/Alias	FACT; T160; FACT80; SSRP1
Species	Human
Gene ID	6749
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human SSRP1 (NP_003137.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	Q08945
Protein Weight	81kDa
Shipping	Ice bag