

Anti-Sstr2 Polyclonal Antibody For Wb, If/Icc, Elisa - P30874

Item Number: CM0024433



Introduction

Rabbit polyclonal antibody targeting human SSTR2 (NP_001041.1), validated for WB, IF/ICC, ELISA; cross-reacts with mouse and rat. Ideal for studying somatostatin signaling and pancreatic hormone regulation.

[Learn More](#)

Parameter	Value
Product Name	Somatostatin Receptor 2 (SSTR2) Rabbit pAb
Remarks/Alias	SST2; Somatostatin Receptor 2 (SSTR2)
Species	Human
GeneID (human)	6752
GeneID	6752
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 270-369 of human Somatostatin Receptor 2 (SSTR2) (NP_001041.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P30874
Protein Weight	41kDa
Shipping	Ice bag