

Anti-Sstr5 Monoclonal Antibody For Wb, Ihc-P, Elisa - P35346

Item Number: CM0005649



Introduction

High-quality rabbit monoclonal antibody targeting human SSTR5, validated for WB, IHC-P, ELISA. Cross-reacts with mouse and rat. Ideal for somatostatin receptor research.

[Learn More](#)

Parameter	Value
Product Name	SSTR5 Rabbit mAb
Remarks/Alias	SST5; SS-5-R; SSTR5
Species	Human
GeneID (human)	6755
GeneID	6755
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human SSTR5 (P35346).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P35346
Protein Weight	39kDa
Shipping	Ice bag