


Anti-Ssx2 Polyclonal Antibody For Wb, If/Icc, Elisa - Q16385

Item Number: CM0013009



Introduction

Rabbit polyclonal antibody against human SSX2 (cancer/testis antigen 5.2), validated for WB, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Ideal for oncology and transcriptional regulation research.

[Learn More](#)

Parameter	Value
Product Name	SSX2 Rabbit pAb
Remarks/Alias	SSX; HD21; CT5.2; CT5.2A; HOM-MEL-40; SSX2
Species	Human
GeneID (Human)	6757
GeneID	6757
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-188 of human SSX2 (NP_783629.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16385
Protein Weight	22kDa
Shipping	Ice bag