

Anti-Ssx4 Polyclonal Antibody For Wb, Elisa - O60224

Item Number: CM0028758

Introduction

SSX4 Rabbit pAb validated for WB and ELISA, targeting human and mouse SSX4. Ideal for transcription regulation research and cancer/testis antigen studies. Supplied by CamelBio, your IVD raw material partner.

[Learn More](#)

Parameter	Value
Product Name	SSX4 Rabbit pAb
Remarks/Alias	CT5.4; SSX4
Species	Human
GeneID (Human)	6759
GeneID	6759
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 70-188 of human SSX4 (NP_005627.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	O60224
Protein Weight	22kDa
Shipping	Ice bag