

Anti-St3Gal1 Rabbit Polyclonal Antibody For Wb And Elisa - Q11201

Item Number: CM0014773

Introduction

High-quality rabbit polyclonal antibody against human ST3GAL1, validated for WB and ELISA, with mouse and rat cross-reactivity. Ideal for glycobiology and sialylation research.

[Learn More](#)

Parameter	Specification
Product Name	ST3GAL1 Rabbit pAb
Remarks/Alias	ST3O; SIAT4A; SIATFL; ST3GalA; ST3GalIA; Gal-NAc6S; ST3GalA.1; 1; ST3GAL1
Species	Human
GeneID	6482
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 35-135 of human ST3GAL1 (NP_003024.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q11201
Protein Weight	39 kDa
Shipping	Ice bag