

Anti-St3Gal4 Polyclonal Antibody For Wb, Elisa - Q11206

Item Number: CM0026033

Introduction

High-quality ST3GAL4 Rabbit Polyclonal Antibody validated for WB and ELISA. Detects human and mouse ST3GAL4, a key sialyltransferase in glycoprotein/glycolipid biosynthesis and hemostasis, with 38 kDa molecular weight.

[Learn More](#)

Parameter	Value
Product Name	ST3GAL4 Rabbit pAb
Remarks/Alias	STZ; SAT3; ST-4; CGS23; SIAT4; NANTA3; SIAT4C; ST3GalIV; ST3GalA.2; gal-NAc6S; ST3GAL4
Species	Human
GeneID (Human)	6484
GeneID	6484
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 27-267 of human ST3GAL4 (NP_006269.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q11206
Protein Weight	38kDa
Shipping	Ice bag