

Anti-St3Gal4 Rabbit Polyclonal Antibody (Wb, Elisa) - Q11206

Item Number: CM0019578

Introduction

High-quality rabbit polyclonal antibody targeting human ST3GAL4 (Q11206). Validated for WB and ELISA, with mouse cross-reactivity. Ideal for studying sialyltransferase function, hemostasis, and leukocyte adhesion.

[Learn More](#)

Parameter	Value
Product Name	ST3GAL4 Rabbit pAb
Remarks/Alias	ST3GAL4; CGS23; NANTA3; SAT3; SIAT4; SIAT4C; ST-4; ST3GalA.2; ST3GalIV; STZ; gal-NAc6S; ST3 beta-galactoside alpha-2; 3-sialyltransferase 4
Species	Human
GeneID (human)	6484
GeneID	6484
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 234-333 of human ST3GAL4(NP_001241686.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q11206
Protein Weight	36kDa/37kDa/38kDa
Shipping	Ice bag