

Anti-St3Gal5 Polyclonal Antibody For Wb, Elisa - Q9Unp4

Item Number: CM0028893

Introduction

Rabbit polyclonal antibody targeting human ST3GAL5 (GM3 synthase), confirmed for WB and ELISA with reactivity across human, mouse, and rat. A valuable tool for ganglioside biosynthesis and neurobiology studies.

[Learn More](#)

Parameter	Value
Product Name	ST3GAL5 Rabbit pAb
Remarks/Alias	SATI; SIAT9; SPDRS; ST3GalV; SIATGM3S; ST3Gal V; ST3GAL5
Species	Human
Gene ID (Human)	8869
Gene ID	8869
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 80-200 of human ST3GAL5 (NP_003887.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UNP4
Protein Weight	48kDa
Shipping	Ice bag