

Anti-St6Gal1 Polyclonal Antibody For Wb, Ihc-P, Elisa - P15907

Item Number: CM0029204



Introduction

Rabbit polyclonal antibody against ST6GAL1 (alpha-2,6-sialyltransferase 1), validated for WB, IHC-P, and ELISA. Detects human, mouse, and rat ST6GAL1, a Golgi enzyme involved in B cell antigen formation and IgG sialylation.

[Learn More](#)

Parameter	Value
Product Name	ST6GAL1 Rabbit pAb
Remarks/Alias	ST6N; CDw75; SIAT1; ST6GalI; ST6GAL1
Species	Human
GeneID (human)	6480
GeneID	6480
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 27-200 of human ST6GAL1 (NP_775323.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P15907
Protein Weight	47kDa
Shipping	Ice bag