

Anti-ST6Gal1/Cd75 Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Fc, Elisa - P15907

Item Number: CM0006422



Introduction

Rabbit monoclonal antibody against ST6GAL1/CD75, validated for WB, IHC-P, IF/ICC, FC, and ELISA. Suitable for studying sialyltransferase function and glycosylation in human samples.

[Learn More](#)

Parameter	Value
Product Name	ST6GAL1/CD75 Rabbit mAb
Remarks/Alias	ST6N; CDw75; SIAT1; ST6Gall; ST6GAL1/CD75
Species	Human
Gene ID (Human)	6480
Gene ID	6480
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 44-406 of human ST6GAL1/CD75 (NP_001340845.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, FC, ELISA
Cross Reactivity	Human
SWISS	P15907
Protein Weight	20kDa/47kDa
Shipping	Ice bag