

Anti-St8Sia1 Polyclonal Antibody For Wb, Elisa - Q92185

Item Number: CM0030986

Introduction

Rabbit polyclonal antibody targeting human ST8SIA1 (GD3 synthase). Detects endogenous ST8SIA1 in WB and ELISA. Cross-reacts with mouse. Ideal for ganglioside metabolism and cancer research.

[Learn More](#)

Parameter	Value
Product Name	ST8SIA1 Rabbit pAb
Remarks/Alias	GD3S; SIAT8; SIAT8A; SIAT8-A; ST8Sial; ST8SIA1
Species	Human
Gene ID (Human)	6489
Gene ID	6489
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 49-270 of human ST8SIA1 (NP_003025.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q92185
Protein Weight	41kDa
Shipping	Ice bag