

Anti-St8Sia2 Polyclonal Antibody For Wb And Elisa - Q92186

Item Number: CM0025087

Introduction

Rabbit polyclonal antibody targeting human ST8SIA2 (alpha-2,8-sialyltransferase 8B), suitable for Western blot and ELISA. Cross-reacts with mouse and rat. Ideal for studying polysialic acid synthesis and neuronal function research.

[Learn More](#)

Parameter	Value
Product Name	ST8SIA2 Rabbit pAb
Remarks/Alias	STX; SIAT8B; SIAT8-B; HsT19690; ST8Siall; ST8SIA-II; ST8SIA2
Species	Human
Gene ID (Human)	8128
Gene ID	8128
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 24-160 of human ST8SIA2 (NP_006002.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q92186
Protein Weight	42kDa
Shipping	Ice bag